



AMI | ITALIA
MEDICAL DEVICES



INFRARED THERMOMETER YK-IRT2



AMI

ITALIA

MEDICAL DEVICES

CONTACTLESS MEASUREMENT
PRECISE E RELIABLE



HANDY, LIGHTWEIGHT
AND EASY TO USE

1 SECOND
RESPONSE TIME!





ACOUSTIC & COLOR INDICATORS



BACKLIT DISPLAY
WITH 3 COLORS
FOR 3
DIFFERENT
TEMPERATURE
RANGES



AMI | **ITALIA**
MEDICAL DEVICES

**LIQUID TEMPERATURE
MEASUREMENT**



SUITABLE FOR BABY

**ENVIROMENT TEMPERATURE
MEASUREMENT**



Product Display



**HANDY
AND EASY TO USE**





PRODUCT INFORMATION

Non-contact digital infrared thermometer with high precision and reliability. Multifunction, handy and easy to use; allows you to detect the temperature in 1 second and to monitor the progress of the fever, by comparing 34 detailed memories with date and time.

FEATURES

- Multiple mode: Body temperature / Surface temperature / Liquid temperature (milk, water)
- 34 detailed memories with date and time, for a simple comparative analysis
- Response time: 1 second
- Measuring distance: up to 5-15 cm, contactless
- Fever alarm (visual and audible)
- Selectable C° / F° mode
- Backlit display with three colors for three different temperature ranges detected
- Parameter setting editable
- Handy and easy to use, with a single multifunction button

USE

- Suitable for use with children
- Suitable for use in a clinical and family environment
- No contraindications

TECHNICAL DATA

- Body temperature range: 34 ~ 42.9 ° C
- Surface temperature range: 0 ~ 100.0 ° C
- Measuring distance: 3-5 cm
- Measurement time: 1 second
- Accuracy: (+/-) 0.2 °
- Power supply: 2 AAA 1.5 V batteries (not included)
- Measurements: 165 x 64 x 66mm
- Dimensions: weight 120 g (without batteries) - measures 15 x 7.5x 5 cm
- Storage and transportation: - 20 - + 45 °C