



AMI | ITALIA  
MEDICAL DEVICES



*MONITORING*

**A.M.I. ITALIA S.r.l.**

Official Distributor

Via Cupa Reginella, 15/A - 80010 QUARTO, NAPLES (ITA)

Tel. +39 081 8060574 - Fax +39 081 8764769

[www.amiitalia.com](http://www.amiitalia.com) - [info@amiitalia.com](mailto:info@amiitalia.com)



## Multi-parameter portable Monitor Model: YK-8000C

### MAIN FEATURES

- ☆ 6 Standard parameters: ECG, RESP, NIBP, SPO2, 2TEMP, PR / HR. 5 Leads: RA, LA, LL, RL, V.
- ☆ 12.1 "high resolution colour LCD display, wide viewing angle
- ☆ Real-time S-T segment analysis, pacemaker detection
- ☆ Indication of handling and malfunction
- ☆ Patient information management function
- ☆ Phase display of multi-conductor ECG waveforms
- ☆ 3-level audio / video alarms
- ☆ Efficient resistance to interference
- ☆ Playback / navigation function for recording waveform and monitor data
- ☆ Built-in rechargeable lithium battery, 2h of operation after full charge

### DIMENSIONS AND CONFIGURATIONS

Standard configuration: Spo2, PR, NIBP, ECG, RESP, TEMP

Optional configuration: Nellcor Spo2, 2-IBP, ETCO2, Touch screen, integrated thermal printer

Dimensions: 30.5 cm x 16.2 cm x 29 cm Weight: 3.3 kg

### ALIMENTAZIONE:

- AC 100-240v AC, 50/60 Hz, maximum input power 70VA
- Battery: 11.1V 2200mAh rechargeable lithium battery

### CONTENT

- 1 YK-8000C device
- 1 lithium battery
- 1 Power cable
- 1 Ground wire
- 1 Probe for Spo2, PR
- 1 Blood pressure measuring components (for NIBP)
- 1 ECG cable (for ECG, RESP)
- 1 Temperature detection probe

**A.M.I. ITALIA S.r.l.**

Official Distributor

Via Cupa Reginella, 15/A - 80010 QUARTO, NAPLES (ITA)

Tel. +39 081 8060574 - Fax +39 081 8764769

[www.amiitalia.com](http://www.amiitalia.com) - [info@amiitalia.com](mailto:info@amiitalia.com)

## Multi-parameter portable Monitor Model: YK-810



### MAIN FEATURES

- ☆ Standard parameters: SPO2, NIBP, PR
- ☆ Options: EtCO2 monitoring, SPO2 neonatal sensor
- ☆ 4.3 " color LCD display, wide viewing angle, high brightness
- ☆ Compact design, comfortably portable
- ☆ Intuitive simple display interface
- ☆ Large capacity internal rechargeable battery
- ☆ Playback and navigation function for long-term data recording.

### DIMENSIONS AND CONFIGURATIONS

Standard configuration: SpO2, NIBP, ECG

Optional configuration: EtCO2, neonatal sensor SPO2, integrated thermal printer

Dimension: 25.5 cm x 15.2 cm x 10, cm Weight: 1.5 kg

### POWER SUPPLY:

-Battery: 11.1V 2000mAh rechargeable lithium

### CONTENT

- 1 YK-810 device
- 1 Lithium battery
- 1 Power cable
- 1 SpO2 probe
- 1 Blood pressure measuring components (for NIBP)
- 1 User manual



## INFORMAZIONI SUL PRODOTTO

Monitor multiparametrico

Modello: E-8

### CARATTERISTICHE PRINCIPALI

- ☆ Parametri standard: ECG, NIBP, SPO2
- ☆ display LCD a colori da 8", ampio angolo di visione, alta luminosità, risoluzione 800\*600
- ☆ Interfaccia di visualizzazione semplice e intuitiva
- ☆ Funzione di riproduzione e navigazione per la registrazione della forma d'onda e dei dati del monitor
- ☆ Doppio allarme automatico con segnali acustici e visivi
- ☆ Visualizzazione in fase delle forme d'onda ECG a più conduttori.
- ☆ Allarmi audio/video a 3 livelli.
- ☆ Efficiente resistenza alle interferenze
- ☆ Efficiente resistenza alle interferenze
- ☆ Batteria al litio ricaricabile incorporata, 2h di funzionamento dopo la carica completa

### DIMENSIONI E CONFIGURAZIONI

Configurazione standard: Spo2, NIBP, ECG

#### ALIMENTAZIONE:

- AC: 100-240V,50Hz/60Hz
- Batteria: 11.1V 4000mAh al litio ricaricabile

### CONTENUTO

- 1 dispositivo E-8
- 1 Batteria al litio
- 1 Cavo di alimentazione
- 1 Cavo di messa a terra
- 1 Sonda per Spo2
- 1 Componenti misurazione Pressione sanguigna (per NIBP)
- 1 Cavo ECG (per ECG, RESP)
- 1 Manuale utente

**A.M.I. ITALIA S.r.l.**

Distributore Ufficiale Italia

[www.amiitalia.com](http://www.amiitalia.com) - [info@amiitalia.com](mailto:info@amiitalia.com)